




Custom accessories to meet your animal research needs

OptiMICE Wire Bar Cage Top Assembly

C79295



Animal Care Systems developed this new accessory specifically for one of our caging customers. The Wire Bar Lid gives OptiMICE users the freedom to house animals in an OptiMICE cage without ventilation from the rack. This allows temporary bench-top animal housing after surgery, or any situation requiring housing animals off the OptiMICE rack. The Wire Bar Lid temporarily replaces the standard OptiMICE cage top, and provides breathable cage lid security when an occupied cage is on the workbench or in an area without active ventilation. Made of electropolished stainless steel, it snaps onto the OptiMICE cage identically to the standard cage top.

Animal Care Systems, Inc., 7086 S. Revere Parkway, Suite 100, Centennial, CO 80112

Phone: 1.888.827.3861 --- Local: 720.283.0177 Fax: 720.283.0179

Contact us: info@animalcaresystems.com

Visit us on the Web: www.animalcaresystems.com