



**Custom accessories to meet your animal research needs**

## **OptiMICE Sac n Pac Feeder**

**C79330H** – for Lab Products Hydropac®  
**C79330S** Sac n Pac Feeder—for Edstrom Sipper Sack®



Feeder is constructed of all stainless steel and fits in the OptiMICE cage in place of the standard feeder and water bottle. It accommodates either the Edstrom Sipper Sack® or the Lab Products Hydropac®.



## **OptiMICE Sac n Pac Solo C79331**

Replaces water bottle in the Sac n Pac.  
Positioned next to the ACS standard feeder.

Animal Care Systems, Inc., 7086 S. Revere Parkway, Suite 100, Centennial, CO 80112  
Phone: 1.888.827.3861 --- Local: 720.283.0177 Fax: 720.283.0179  
Contact us: [info@animalcaresystems.com](mailto:info@animalcaresystems.com)  
Visit us on the Web: [www.animalcaresystems.com](http://www.animalcaresystems.com)

Animal Care Systems, Inc., 7086 S. Revere Parkway, Suite 100, Centennial, CO 80112  
Phone: 1.888.827.3861 --- Local: 720.283.0177 Fax: 720.283.0179  
Contact us: [info@animalcaresystems.com](mailto:info@animalcaresystems.com)  
Visit us on the Web: [www.animalcaresystems.com](http://www.animalcaresystems.com)