



**Custom accessories to meet your animal research needs**

## **Opti Rack Status Monitor**

**C79351**



The C79351 Rack Status Monitor provides continuous rack airflow monitoring on all OptiMICE, ErgoMICE and OptiRAT racks. Simple hookup and battery power ensure quick and easy installation on any Opti rack. Shipped complete with display hanger and batteries.

Animal Care Systems, Inc., 7086 S. Revere Parkway, Suite 100, Centennial, CO 80112  
Phone: 1.888.827.3861 --- Local: 720.283.0177 Fax: 720.283.0179  
Contact us: [info@animalcaresystems.com](mailto:info@animalcaresystems.com)  
Visit us on the Web: [www.animalcaresystems.com](http://www.animalcaresystems.com)